

## INNOCONSUMABLE™ General Laboratory Consumables Series

### General Description

The consumables provided by Velsson are compatible with all Velsson automated equipment and are single-use products. These consumables are sterilized by gamma irradiation and tested to be pyrogen-free. They are also compatible with a range of third-party laboratory equipment.

### Product Features

- Free of human DNA, RNase, and DNase.
- Shelf life of 1 year.
- Non-cytotoxic and pyrogen-free.
- Free of BSE/TSE risk.
- Compliant with the following international standards:
  - United States Pharmacopeia (USP VI)
  - European Pharmacopoeia (EP)
  - ISO Standards (ISO 10993-1, 13485:2016)
- Sterilized by gamma irradiation:
  - Minimum gamma dose: 15.0 kGy
  - Maximum gamma dose: 25.0 kGy
  - Endotoxin level: < 0.1 EU/mL

### Product Storage

Proper storage conditions are essential to maintaining the integrity and performance of Velsson products.

- Store at room temperature.
- Avoid UV or direct sunlight exposure.

### Warning

Velsson products are designed to achieve optimal performance under specific conditions. Do not use the product if:

- The packaging has been damaged.

### Consumables Comparison Table

The following table provides a general comparison of Velsson's various products, highlighting key features and specifications to help you select the right consumable for your application.

## Tubes

Item	Product No.	Per Case
50 mL Centrifuge Tube	VEL-CS-002001	500 pcs/case
50 mL Centrifuge Tube – 50 mL Auto-Lock Culture Tube	VEL-CS-002002	500 pcs/case
50 mL Culture Tube (500 pcs/case)	VEL-CS-002003	500 pcs/case
50 mL Culture Tube	VEL-CS-002003	500 pcs/case
50 mL Auto-Lock Culture Tube	VEL-CS-002004	500 pcs/case

## Pipette Tips

Item	Product No.	Per Case
10 µL Tips and Tip Box	VEL-CS-002014	10 boxes/rack, 5 racks/case
50 µL Tips and Tip Box	VEL-CS-002015	10 boxes/rack, 5 racks/case
300 µL Tips and Tip Box	VEL-CS-002016	10 boxes/rack, 5 racks/case
300 µL Extended Tips and Tip Box	VEL-CS-002034	10 boxes/rack, 5 racks/case
1000 µL Tips and Tip Box	VEL-CS-002017	10 boxes/rack, 5 racks/case

## Reagent Reservoirs / Shake Flasks

Item	Product No.	Description	Per Case
Reagent Reservoir	VEL-CS-002005	Square well, V-bottom, 300 mL capacity	50 pcs/case